

FLJ34048-hypothetical transcript Gene

Reference: PTXBC113049

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical transcript
GENE ID	285987
ORIGIN SPECIES	Human
ACCESSIONS	BC113049
VECTOR	pDONR223
NCBI SYMBOL	FLJ34048
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggactggggaagacattaaagtcggagcagcccaaggctcggatctagcctttctccctccggtctccagacccacgcaggttaaggcagccaagcagagacagcggagccccgggtcgtggttcaggcctggtcaacaggtaggaggccagactggggccaggcacag
aggaggcaggccccccgacccggcgtccagggtccccgggtccgggtccgggactggaggaggtcggcgagccccccctccagaggcttgagattgggtggagacagctgaagatcctctcgtctggaccccgactttcgatctgcttag