

FLJ38359-hypothetical protein FLJ38359 Gene

Reference: PTXBC113052

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC151009		

For research use only.

GENE DESCRIPTION	hypothetical protein FLJ38359
GENE ID	151009
ORIGIN SPECIES	Human
ACCESSIONS	BC113052
VECTOR	pDONR223
NCBI SYMBOL	FLJ38359
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacacagtcacgtcgacagaagtcccaccagctctgcccattgaggggtctcaagcagagagcgtgagaagttctgcagatgtatgtggcgacagcctctacagccgagctccgactctcttcaggcaccaggcctggccgtggccgctctcaccgactccaaggcagc
actccggcgctctctttattttatcacagttttatgagatataatttacatgctatccagcccatctattcaagtgtacatgcaagggtttgtatataatcacagagttgtgtatccatcaccatcgttttagaacattcactcactccagaagaactactacccatcagcagccactgcccatttc
cctga