

KRTAP3-1-keratin associated protein 3-1 Gene

Reference: PTXBC113077

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP3.1; KRTAP3.1; keratin-associated protein 3-1; high sulfur keratin-associated protein 3.1; keratin-associated protein 3.1; keratin associated protein 3-1		

For research use only.

GENE DESCRIPTION	keratin associated protein 3-1
GENE ID	83896
ORIGIN SPECIES	Human
ACCESSIONS	BC113077
VECTOR	pDONR223
NCBI SYMBOL	KRTAP3-1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtattgtgtgtctctcgcctcgcagcgtccccaccggccctgccaccaccttctgctcatttgataaaagctgccgctgtggagctgtgctaccagcaccctgccacatgagatcagcctcttcagcccatctgtgtgacacctgccccccaccctgtgcaagcctgatactatgt
gcccaactgtgtgtgtctcaacaactgtcacccgactccggactgagtgaggatcaacctgaccacctatgttcagcctggctgtgagagtcctgtgagccccgctgttaa