

SYPL2-synaptophysin-like 2 Gene

Reference: PTXBC113103

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MG29; synaptophysin-like protein 2; mitsugumin 29; synaptophysin like 2		

For research use only.

GENE DESCRIPTION	synaptophysin-like 2
GENE ID	284612
ORIGIN SPECIES	Human
ACCESSIONS	BC113103
VECTOR	pDONR223
NCBI SYMBOL	SYPL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctctgcgatgaagagtcacgtccaagaccatgcacctcatgggggacttctctgcacccgccgagttctctgtgaccttggcatctttctctctataccatggctgccctgtatctacgtgcgttccacaacctctacacagagaacaaacgttcccgctggtggacttctgtgtgactgtctcttcccttctctgtgctggtgagctgcagctgcctggggcaaggccctgaccgatgtcaaggggccacacgacctccagctgacagcagcatgtcagttgccatggagaggagcagtggtgagtgccggggccacgccctctatggcctggccaacatctccgtggtgagacctgtggccactgcaggagcagcaccagccctgctgcccaggcctgtcccagctag