



Reference: PTXBC113370

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NCRNA00153; C20orf19; HT013; Kizuna; PLK1S1; RP69; centrosomal protein kizuna; polo-like kinase 1 substrate 1; kizuna centrosomal protein		

For research use only.

GENE DESCRIPTION	non-protein coding RNA 153
GENE ID	55857
ORIGIN SPECIES	Human
ACCESSIONS	BC113370
VECTOR	pDONR223
NCBI SYMBOL	NCRNA00153
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr