

FLJ32063-hypothetical protein LOC150538 Gene

Reference: PTXBC113374

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC150538
GENE ID	150538
ORIGIN SPECIES	Human
ACCESSIONS	BC113374
VECTOR	pDONR223
NCBI SYMBOL	FLJ32063
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtatggctctataccacacagcttattctacaagaagtccccgaggaggattgggttaagctttgcaaatgtgcttccaggtaatgcgttcattattctgctccgacttaccaccacaccagtggtaccggagcagcaccacttggcagaactgatgactgcttggccca
gcggagtgccgattgcgctaacacgcgcaggggaattgcacccttccggagcctccgcacgtgctgcccctcaaaagagctggcgaccggctcacgtgtaagcaacctccactttgaaactaatcgacccgggtctttaccccaaggactttgctgcggacgtgctctgaccc
aagacgcgggagagaagtcccaaggctacagccaggggctgggggactcctctgctcaccttagtccttgatttgaaggcccaattaa