



Reference: PTXBC113381

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	histone demethylase UTX; bA386N14.2 (ubiquitously transcribed X chromosome tetratricopeptide repeat protein (UTX)); UTX; KABUK2; bA386N14.2; lysine-specific demethylase 6A; lysine (K)-specific demethylase 6A; ubiquitously transcribed tetratricopeptide rep		

For research use only.

GENE DESCRIPTION	ubiquitously transcribed tetratricopeptide repeat, X chromosome
GENE ID	7403
ORIGIN SPECIES	Human
ACCESSIONS	BC113381
VECTOR	pDONR223
NCBI SYMBOL	UTX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]