



Reference: PTXBC113393

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HTX; HTX1; VACTERLX; ZNF203; zinc finger protein ZIC 3; Zic family member 3 (odd-paired homolog, Drosophila); heterotaxy 1; zinc finger protein 203; zinc finger protein of the cerebellum 3; Zic family member 3		

For research use only.

GENE DESCRIPTION	Zic family member 3 (odd-paired homolog, Drosophila)
GENE ID	7547
ORIGIN SPECIES	Human
ACCESSIONS	BC113393
VECTOR	pDONR223
NCBI SYMBOL	ZIC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1