



## Reference: PTXBC113402

|                |                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | Fz2; fz-2; fzE2; hFz2; frizzled-2; frizzled 2, seven transmembrane spanning receptor; frizzled family receptor 2; frizzled homolog 2; frizzled class receptor 2 |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | frizzled homolog 2 (Drosophila)                                                                   |
| <b>GENE ID</b>          | 2535                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC113402                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | FZD2                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

[illegible]

More info: Please contact us for any question: [contact@proteogenix.fr](mailto:contact@proteogenix.fr)