



Reference: PTXBC113406

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ICE(rel)III; ICEREL-III; ICH-3; caspase-5; CASP-5; ICE(rel)-III; TY protease; caspase 5, apoptosis-related cysteine peptidase; caspase 5, apoptosis-related cysteine protease; protease ICH-3; protease TY; caspase 5		

For research use only.

GENE DESCRIPTION	caspace 5, apoptosis-related cysteine peptidase
GENE ID	838
ORIGIN SPECIES	Human
ACCESSIONS	BC113406
VECTOR	pDONR223
NCBI SYMBOL	CASP5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr