



Reference: PTXBC113423

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATS; GLUT10; solute carrier family 2, facilitated glucose transporter member 10; GLUT-10; glucose transporter type 10; solute carrier family 2 (facilitated glucose transporter), member 10; solute carrier family 2 member 10		

For research use only.

GENE DESCRIPTION	solute carrier family 2 (facilitated glucose transporter), member 10
GENE ID	81031
ORIGIN SPECIES	Human
ACCESSIONS	BC113423
VECTOR	pDONR223
NCBI SYMBOL	SLC2A10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr