



Reference: PTXBC113427

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZABC1; zinc finger protein 217		
For research use only.			

GENE DESCRIPTION	zinc finger protein 217
GENE ID	7764
ORIGIN SPECIES	Human
ACCESSIONS	BC113427
VECTOR	pDONR223
NCBI SYMBOL	ZNF217
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr