

**UBE2J2-ubiquitin-conjugating enzyme E2, J2 (UBC6 homolog, yeast) Gene**

Reference: PTXBC113430

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NCUBE-2; NCUBE2; PRO2121; ubiquitin-conjugating enzyme E2 J2; E2 ubiquitin-conjugating enzyme J2; non-canonical ubiquitin-conjugating enzyme 2; ubiquitin conjugating enzyme 6; ubiquitin-conjugating enzyme E2, J2 (UBC6 homolog, yeast); yeast UBC6 homolog;		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	ubiquitin-conjugating enzyme E2, J2 (UBC6 homolog, yeast)
<b>GENE ID</b>	118424
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC113430
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	UBE2J2
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atgacgagcaccagcagtaagaggctccgaccacggcaaccagaggctgaagcaggactacctcgcattaagaaagaccggctgcttacatctgtgcccagcccctccctcgaatattctcagtgaggcactatgtctccgagggccagagatgacccttatgaaggtgctat  
tatcatggaaaactaattttccagagaattcttcaaacctccagtatctatgatactccaacggagggttaagtgaacaccaggctgtgtcttctatcacggattccaccggacacgtggaacccggcctgtgtctcaccatcctgactgggctcctgagctcatggt  
ggagaaggcccccacctgggcagtagagacgtcggactcacgaaagacaactggcagtgacagtttagcafttaattgaaagataaagcttttgaattatctcgaagtcgtgagagagattaacaaaaacagaaacacaagacgaactcagtagcagaccagact  
ctccccttgcagacgtggttccagacgggagacgcacttgtccagacgggattcagctgctcaacggcctgcgccggggccctcccaacctcgcaggctccagcagccaaccggcaccacggactcctgggtggcgcctggcgaactgtttgtatagttgggtttg  
cagcctttgcttacacgggtcaagtactgctgaggagcctcgcagagtgta