



Reference: PTXBC113438

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GPCR135; RLN3R1; RXFPR3; SALPR; relaxin-3 receptor 1; G-protein coupled receptor GPCR135; G-protein coupled receptor SALPR; RLN3 receptor 1; g protein-coupled receptor SALPR; relaxin family peptide receptor 3; somatostatin and angiotensin-like peptide rec		

For research use only.

GENE DESCRIPTION	relaxin/insulin-like family peptide receptor 3
GENE ID	51289
ORIGIN SPECIES	Human
ACCESSIONS	BC113438
VECTOR	pDONR223
NCBI SYMBOL	RXFP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]