



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MIXL1-Mix1 homeobox-like 1 (*Xenopus laevis*) Gene

Reference: PTXBC113441

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	homeobox protein MIXL1; MILD1; MIX; MIXL; MIX1 homeobox-like protein 1; Mix-like homeobox protein 1; homeodomain protein MIX; mix.1 homeobox-like protein; Mix paired-like homeobox		

For research use only.

GENE DESCRIPTION	Mixl1 homeobox-like 1 (<i>Xenopus laevis</i>)
GENE ID	83881
ORIGIN SPECIES	Human
ACCESSIONS	BC113441
VECTOR	pDONR223
NCBI SYMBOL	MIXL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccacagccgagctccgctgctcagtttgcgaaggccgcgcgtttcacgctgacacggggccccacacgcgcgcggggcgctctctgcgcgcccgagccctgcggcagccctgctcctgcgcgcccgccggcccccggccagcagcttctcggcttctctgcgcgg
gaccctggggccggccccgcgcgcccccgcagcttggggctgctctgcgccccccaaaagcgcggcgccgcctgcctgcgcgtcgcagcggccgcaagcgcacgttcttcagcgcgcgaacagctgcagctgctggagctgctcttctccgcgggaccgggtaccgccagatcacttgcgcg
gagcgcctgcgcgcgcacccctgctccccgagctcggagatccagctatgcttgcgttcacagaacagctgcgcgaagtgcagctgagctgggaatctcgaacacttcttgctclaggcgcgcggagattatctcaaacactgctgctctggctgaacacgaatgttgcgcg
ttgagctagatgtgaactctccgcgcgaacaaacgggggttgaggggcggctctgcactctgacatctcccaagctcagaatttgaaccttgcacccctctctcgaacgtgttcgaacgctgacactgcgcgaggaagaacacttcttgccttgccttgcactacttga

More info: Please contact us for any question: contact@proteogenix.fr