

LCE4A-late cornified envelope 4A Gene

Reference: PTXBC113446

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LEP8; SPRL4A; late cornified envelope protein 4A; late envelope protein 8; small proline rich-like (epidermal differentiation complex) 4A; small proline-rich-like epidermal differentiation complex protein 4A; late cornified envelope 4A		

For research use only.

GENE DESCRIPTION	late cornified envelope 4A
GENE ID	199834
ORIGIN SPECIES	Human
ACCESSIONS	BC113446
VECTOR	pDONR223
NCBI SYMBOL	LCE4A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcctgccagcagaaccaacagcagtgccagccctcccaagtgtctatcccaagtatcccccaaatgtccctcaaatgtgcacatcatgccacctccaatctcttctgtgtggtccagctctgggggctgtggtgtcgcagctctgaggaggtggtgctgcctgagc
caccacagacaccatagggtccactgccacagacccaagagctccaattgctatggcagtgccagtgccagcagctctgggggtctgtgctgtctctggaggggctgtgtctga