



Reference: PTXBC113449

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM9A; TEX39A; testis expressed 39A; family with sequence similarity 9 member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 9, member A
GENE ID	171482
ORIGIN SPECIES	Human
ACCESSIONS	BC113449
VECTOR	pDONR223
NCBI SYMBOL	FAM9A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]