

PTPRB-protein tyrosine phosphatase, receptor type, B Gene

Reference: PTXBC113463

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HPTP-BETA; HPTPB; PTPB; R-PTP-BETA; VEPTP; receptor-type tyrosine-protein phosphatase beta; VE-PTP; protein tyrosine phosphatase, receptor type, beta polypeptide; vascular endothelial protein tyrosine phosphatase; protein tyrosine phosphatase, receptor ty		

For research use only.

GENE DESCRIPTION	protein tyrosine phosphatase, receptor type, B
GENE ID	5787
ORIGIN SPECIES	Human
ACCESSIONS	BC113463
VECTOR	pDONR223
NCBI SYMBOL	PTPRB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr