



Reference: PTXBC113469

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WHSC1L1; KMT3F; KMT3G; WHISTLE; pp14328; histone-lysine N-methyltransferase NSD3; WHSC1-like 1 isoform 9 with methyltransferase activity to lysine; Wolf-Hirschhorn syndrome candidate 1-like 1; nuclear SET domain-containing protein 3; protein whistle; nucl		

For research use only.

GENE DESCRIPTION	Wolf-Hirschhorn syndrome candidate 1-like 1
GENE ID	54904
ORIGIN SPECIES	Human
ACCESSIONS	BC113469
VECTOR	pDONR223
NCBI SYMBOL	WHSC1L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]