



Reference: PTXBC113471

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PVRL1; CD111; CLPED1; ED4; HlgR; HV1S; HVEC; OFC7; PRR; PRR1; PVRR; PVRR1; SK-12; nectin-1; ectodermal dysplasia 4 (Margarita Island type); herpes simplex virus type 1 sensitivity; herpes virus entry mediator C; herpesvirus Ig-like receptor; nectin 1; pol		

GENE DESCRIPTION	poliovirus receptor-related 1 (herpesvirus entry mediator C)
GENE ID	5818
ORIGIN SPECIES	Human
ACCESSIONS	BC113471
VECTOR	pDONR223
NCBI SYMBOL	PVRL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]