

DAND5-DAN domain family, member 5 Gene

Reference: PTXBC113476

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CER2; CERL2; CKTSF1B3; COCO; CRL2; DANTE; GREM3; SP1; DAN domain family member 5; DAN domain family member 5, BMP antagonist; DAN domain family, member 5; cerberus 2; cerberus-like 2; cerberus-like protein 2; cerl-2; cysteine knot superfamily 1, BMP antag		

For research use only.

GENE DESCRIPTION	DAN domain family, member 5
GENE ID	199699
ORIGIN SPECIES	Human
ACCESSIONS	BC113476
VECTOR	pDONR223
NCBI SYMBOL	DAND5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccttggccagctatcactcttctgtgcctgttagcggggccctgctacaggtcagggaggcctgaacccagttctctgacctcagtcctgggtgcagccaatcagacctgggtctgggccaggggccctgccccactggtgccagcttctgacctgggagctgga
aggccttcttgggctgcagaaagccaggcagctggggatgggcaggtgcagcgtgggcaagacgaggtggctgtgtgactctgccgtgaacctcagggaagtgcaggggatgttaaggctgtgccccctgttcagggtttctccggccccgctgtcagccatacgcctc
cgaatcatctgtgcttggctcctctctacatccctggctcgaccccaaccactagtcctgtgcaacagctgtatgcctgtcgaagcgttgggcacccgtgtctctgtgtgtctcactggcagctcagctcccgctgcaggggtgaagatatccacctgctgacgagg
gtgtcactgcagccaaaagcatga