

INDOL1-indoleamine-pyrrole 2,3 dioxygenase-like 1 Gene

Reference: PTXBC113496

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	INDOL1; indoleamine 2,3-dioxygenase 2; IDO-2; indoleamine 2,3-dioxygenase-like 1 protein; indoleamine 2,3-dioxygenase-like protein 1; indoleamine-pyrrole 2,3 dioxygenase-like 1; indoleamine-pyrrole 2,3-dioxygenase-like protein 1		

For research use only.

GENE DESCRIPTION	indoleamine-pyrrole 2,3 dioxygenase-like 1
GENE ID	169355
ORIGIN SPECIES	Human
ACCESSIONS	BC113496
VECTOR	pDONR223
NCBI SYMBOL	INDOL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgcatcttcattatgatactcaacaaaataatggagcccccagaccgaatgtgaagacagcagtgccattgtcttggaaagctatcacatatctgaagatgatggcttcttctccagattctctgaaagaactccagatcattataggccttggatggaaattgccaacaaacttc
ctcaattgattgatgctcaccagctcaagctcatgtggacaagatgccccctgctgagctgccagttctgaagggtcaccgggagcagcgccctggccacctggtcctgagcttctcaccatgggttatgtctggcaggaaggagagcgagcctgcagaggtctgccaagggaatct
tgcccttcatttgcgaagctccaggaaactggggctccctctatctggtccactcagacttgggtgctgacgaactggacaaaaaagatccagacggagacggagctccctctgtctccaggctggagtgcatggtgcatga