

RAB40A-RAB40A, member RAS oncogene family Gene

Reference: PTXBC113501

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB40A, member RAS oncogene family; RAR2; RAR2A; ras-related protein Rab-40A; SOCS box-containing protein RAR2A; protein Rar-2		

For research use only.

GENE DESCRIPTION	RAB40A, member RAS oncogene family
GENE ID	142684
ORIGIN SPECIES	Human
ACCESSIONS	BC113501
VECTOR	pDONR223
NCBI SYMBOL	RAB40A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcggccggcagccccaccagccctatgactctctgctcaagtctctgctggggcgacaggagctaggcaagagtgagatcctggagagcctgcaggatggtgcagctgagtcctccgtacagccatctcgggggatcgactacaagacgaccatcctgctggac
ggccagcgggtgaagctgaagctctgggatacgtcggggcagggaagattttgacatattccgctcctactctcgtggtgcacaaggagtgatcctgggtctacgacattgcaaacgctggtcttgcagggtatggatcgtggaataagaagattgaggaacatgccctggtgctcct
aaaatcctggtggggaatgcctacatctggcattcaagaggcaggtgcccaggagcagggcccgctacgcccagcgcctgggtgtagctcttctgaggtcagccctctgcaatttcaacatcatagagctcttcacggagctggccaggatgctgctgctggcagaggatg
aattggctcgggagccgagcaaggtactgagcttgcagacctctgctgccgacacatcgtgctcgtgcacacctgtgcatctggtggacaagctcccgtccccagtlaccttaagaagccacctcaagtcctctccatggctaaggcctgaatgccagatgctgagggcctccta
ctccctcaccacagctccactcacaagagcagcctctgcaagtgaggatcgtgctccaccgagagcccccaccccaaaactgcaccagaacagctgcaaaattcttaa