

RETNLB-resistin like beta Gene

Reference: PTXBC113502

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FIZZ1; FIZZ2; HXCP2; RELM-beta; RELMb; RELMbeta; XCP2; resistin-like beta; C/EBP-epsilon regulated myeloid-specific secreted cysteine-rich protein precursor 2; colon and small intestine-specific cysteine-rich protein; colon carcinoma-related gene protein;		

For research use only.

GENE DESCRIPTION	resistin like beta
GENE ID	84666
ORIGIN SPECIES	Human
ACCESSIONS	BC113502
VECTOR	pDONR223
NCBI SYMBOL	RETNLB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggcgctctctgctctctctatcctaataccccctctccagctgatcaacccggggagtagctcagtggttcttagactccgttatggataagaagatcaaggatgttctcaacagctagagtagacagtcctctctataagcaagaagctctctgtgtctagtgtaaaagccaaggca
gaccgtctctctgcccctgctgggatggctgtcactggctgtgcttggctatggctgtggttcgtgggatgttcagctggaaaccacctgccactgccagtcagtggtggactggaccactgcccgtgctgccactgacctga