

BCL2L2-BCL2-like 2 Gene

Reference: PTXBC113522

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BCL-W; BCL2-L-2; BCLW; PPP1R51; bcl-2-like protein 2; apoptosis regulator BCL-W; protein phosphatase 1, regulatory subunit 51; BCL2 like 2		

For research use only.

GENE DESCRIPTION	BCL2-like 2
GENE ID	599
ORIGIN SPECIES	Human
ACCESSIONS	BC113522
VECTOR	pDONR223
NCBI SYMBOL	BCL2L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgacccagcctcgccccagacacacgggctctgggtggcagacttttaggtataagctgaggcagaagggttatgtctgtgagctggccccggggagggccagcagctgacccactgcaccaagccatcggggagctggagatgagttcgagacccgttcggcg
caccttctctgatctggcgctcagctcatgtgacccaggtcagcccaacaacgcttcaccaggtctcgcgatgaacttttcaaggggcccaactggggccgctgtgaccttctttgtcttggggctgactgtgtctgagagtgtaacaaggagatggaaccactgtggg
acaagtgcaggagtgatgggtggcctacctggagacgcggtggtgactggatccacagcagtgggggctggcgagggttcacagctctatagggggacgggcccctggaggagggcgcgctgtctggggaggggaactgggcatcagtgaggacagtgctgacggggccg
tggcactgggggcccctgtaactgtaggggcctttttgtagcaagtga