



Reference: PTXBC113525

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	neuroligin-4, Y-linked; neuroligin 4, Y-linked		
For research use only.			

GENE DESCRIPTION	neuroligin 4, Y-linked
GENE ID	22829
ORIGIN SPECIES	Human
ACCESSIONS	BC113525
VECTOR	pDONR223
NCBI SYMBOL	NLGN4Y
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr