

LOC340094-hypothetical LOC340094 Gene

Reference: PTXBC113532

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC340094
GENE ID	340094
ORIGIN SPECIES	Human
ACCESSIONS	BC113532
VECTOR	pDONR223
NCBI SYMBOL	LOC340094
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacctgcagtgacctgctgctaccagaactgactcgttcaaaagtacggatgaacccataggagatattctgggtgggtaggccagttggaagacagagtgtttacacctgagagaggacttgagaattttactcagtcagttgctcactgctgcgcactcctgcagacgcaccc
cagttggactctgcgctgcatccccagagtacagcgagtttccgcctccgttagagatccaacacatatttttgaagaagaactgatggaagacaggaagcaaaaaggagaggaaagggaaggataggaaggaggggaaggaaagggaaggaaagggaagg
tgaactgaaatcctacaaggctcaaaagtcacaaacatgattgctgttttagagaattgaagatcctcacttcaagaatctgaaatctaccatgatctcacatgtag