

CHMP4A-chromatin modifying protein 4A Gene

Reference: PTXBC113533

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf123; CHMP4; CHMP4B; HSPC134; SHAX2; SNF7; SNF7-1; VPS32-1; VPS32A; charged multivesicular body protein 4a; SNF7 homolog associated with Alix-2; Snf7 homologue associated with Alix 2; chromatin modifying protein 4A; vacuolar protein sorting-associate		

For research use only.

GENE DESCRIPTION	chromatin modifying protein 4A
GENE ID	29082
ORIGIN SPECIES	Human
ACCESSIONS	BC113533
VECTOR	pDONR223
NCBI SYMBOL	CHMP4A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgcggcgccctgaggacggattgggcaagctgtccctgtgtgatgagacatcacctcccaggagcaaggcggaagtctggaggacgtgaggcgccggaggcgagagcgagctcgcatgagtggtctcggcaggctctcgggaagggaaggagaa
aggccaacccctgaagaagcaatacagaactgaaggagacagagaagatactgatcaagaacagggaattttggagcagaagattcaacaggagctacaacagccaagaagtatgggaccaagaataagagagctgccctacaggctttgcggaggagaaaagattcgaac
agcagctggcacaactacgggacattatccacctggagtttcagctgagggccattgagaatgccactaccaatgcagaagtctctgtaccatggagcttctgcccgaagcatgaagaaggcctaccaggacatggacattgacaaggtagatgaactgatgactgacatcacgg
aacaacaggaggtggcccagcagatctcagatccattctcggcctatgggctttggagatgatgtggatgaggatgaactgtggaggagctagaggagctggagcaggaggaattggccaggaggtgttaaatgtggcgacaaggagaagaaccctcagtaaatgcctagt
gtaccttctactcatctgccggcaggccagctcccaaaagtggatgaagatgaagaagcactaaagcagttggctgagtggtatcctga