

KRTAP13-1-keratin associated protein 13-1 Gene

Reference: PTXBC113536

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP13.1; KRTAP13.1; keratin-associated protein 13-1; high sulfur keratin associated protein 13.1; keratin associated protein 13-1		

For research use only.

GENE DESCRIPTION	keratin associated protein 13-1
GENE ID	140258
ORIGIN SPECIES	Human
ACCESSIONS	BC113536
VECTOR	pDONR223
NCBI SYMBOL	KRTAP13-1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggcatgcaaaactcagaatctctcagtgtaactcagctgaactcacatccccatcaacatgtcctacaactgctgcttggaactctctcccgctcctgtgggtggtacctgactaccagcctcctcctgtggctttctaccaccagcaaccagggtctacagcactgacctgtgctc
tcccagcagtgccagctgggttctctctctatagggctgtcagcagacctgctggagcccaccagctgccagacatcctatgtggagtcagccctgccagacctcctgtaccgtcccagaacctcctgtctctgagtcctgccagacaactactctgggtctctaggtttgg
atccagcagctgccgtccctgggctatggatcgaggagctgctactcagtggtggtggtccagtggttcagatccctgggttatggaggtgtggtctccctccctgggctatggcgttggttctgcccgaacctacttggttctaggagctgccagttcttctgctacagacaa
cttggatcaggtcttctactattga