



Reference: PTXBC113547

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GHDP; growth hormone secretagogue receptor type 1; GH-releasing peptide receptor; GHRP; GHS-R; ghrelin receptor		

For research use only.

GENE DESCRIPTION	growth hormone secretagogue receptor
GENE ID	2693
ORIGIN SPECIES	Human
ACCESSIONS	BC113547
VECTOR	pDONR223
NCBI SYMBOL	GHSR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]