



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

GAL3ST3-galactose-3-O-sulfotransferase 3 Gene

Reference: PTXBC113552

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GAL3ST-3; GAL3ST2; galactose-3-O-sulfotransferase 3; beta-galactose-3-O-sulfotransferase 3; gal-beta-1, 3-GalNAc 3'-sulfotransferase 3; galactose 3'-sulfotransferase; galbeta1-3GalNAc 3'-sulfotransferase 3		

For research use only.

GENE DESCRIPTION	galactose-3-O-sulfotransferase 3
GENE ID	89792
ORIGIN SPECIES	Human
ACCESSIONS	BC113552
VECTOR	pDONR223
NCBI SYMBOL	GAL3ST3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr