



Reference: PTXBC113559

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	P4-ATPase flippase complex beta subunit TMEM30B; cell cycle control protein 50B; transmembrane protein 30B		

For research use only.

GENE DESCRIPTION	transmembrane protein 30B
GENE ID	161291
ORIGIN SPECIES	Human
ACCESSIONS	BC113559
VECTOR	pDONR223
NCBI SYMBOL	TMEM30B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1