



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

TAS2R7-taste receptor, type 2, member 7 Gene

Reference: PTXBC113561

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	T2R7; TRB4; taste receptor type 2 member 7; taste receptor, family B, member 4; taste receptor, type 2, member 7; taste 2 receptor member 7		

For research use only.

GENE DESCRIPTION	taste receptor, type 2, member 7
GENE ID	50837
ORIGIN SPECIES	Human
ACCESSIONS	BC113561
VECTOR	pDONR223
NCBI SYMBOL	TAS2R7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr