

C17orf44-chromosome 17 open reading frame 44 Gene

Reference: PTXBC113564

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C17orf44; NCRNA00324; long intergenic non-protein coding RNA 324		
For research use only.			

GENE DESCRIPTION	chromosome 17 open reading frame 44
GENE ID	284029
ORIGIN SPECIES	Human
ACCESSIONS	BC113564
VECTOR	pDONR223
NCBI SYMBOL	C17orf44
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgattctcgcttcgggtgcgagaggtccgggttcaaatccggagagccctcttacctttactgagacaagagtgtctcaaggaatatgtacctgaccagcgtgcaggagcgttgactgaatcatcaatgtacagtaagctctttgtggcaccggttagcagaaagcatctc
aatccccagaccgctggcatgtcagggaagagacgagcagccagacttggcagggtccctcttctgtctcacttctcctgtgtctcaactcggccaccacagcctgtttcgtttgattgacacgcgtcaatctatagtgtggatgacagtgttcgggagagccagga
actgtcaaaacttggttctctgcagctctgaaggcactctgtctcttcggtcactccatctgacgattggtgcctga