

LOC152225-hypothetical LOC152225 Gene

Reference: PTXBC113578

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC152225		

For research use only.

GENE DESCRIPTION	hypothetical LOC152225
GENE ID	152225
ORIGIN SPECIES	Human
ACCESSIONS	BC113578
VECTOR	pDONR223
NCBI SYMBOL	LOC152225
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaataatgtttcaataaggaagacagaatgagttctgacacgatgtgtggagctgtgacaggcagacgaaaaatggagccaaatggcatggaggtgtgtcatcactgctggattttaccttgatctacattcagttatctactcttttcaaaatgctgggcattcattcaagaacagcatat
atgttctgactttgaagtaatggatgaattgtcatgtgcagtatatggaaataagttttattaccttcttccacatttaactcatccatccattcaatag