

IL28A-interleukin 28A (interferon, lambda 2) Gene

Reference: PTXBC113581

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL28A; IL-28A; interferon lambda-2; IFN-lambda-2; cytokine Zcyto20; interleukin 28A (interferon, lambda 2); interleukin-28A; interferon lambda 2		

For research use only.

GENE DESCRIPTION	interleukin 28A (interferon, lambda 2)
GENE ID	282616
ORIGIN SPECIES	Human
ACCESSIONS	BC113581
VECTOR	pDONR223
NCBI SYMBOL	IL28A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaactagacatgactgggactgcacgccagtgtgtgtgatggcgcagtgctgacctgactggagcagttcctgtcgcaggtccacggggctctccggatgcaaggggctgccacatagcccagttcaagtcctgtctccacaggagctgcaggccttaagaggc
caaagatgccttagaagagtcgttctgtgaaggactgcaggtgccactcccgcctcttccccaggacctgggacctgaggcagctgcaggtgagggagcgcgccatggctttggagctgagctggccctgacgtgaaggtcttgaggccaccgtgacactgacccagccctg
gtggacgtctggaccagccccctcacacctgcaccatactctccagttccggcctgtatccagctcagccacggcagggccagggacccggggccgctccacattggctgtaccggctccaggaggcccaaaaaggagttccctggctgcctcagggcctctgtcac
cttcaacctctccgctcctcagcgagacctgaattgtgtgccagtgaggacctgtgtgtctga