

MPZL3-myelin protein zero-like 3 Gene

Reference: PTXBC113586

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	myelin protein zero-like protein 3; myelin protein zero like 3		
For research use only.			

GENE DESCRIPTION	myelin protein zero-like 3
GENE ID	196264
ORIGIN SPECIES	Human
ACCESSIONS	BC113586
VECTOR	pDONR223
NCBI SYMBOL	MPZL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagcagagaggagcagctggaagccgtggctgcgtctcttccctctgctggcgctctgttcttcaggggtgtttatcgtcttcttggagattcgtgcagatgcccatgtccgaggttatgttgagaaaagatcaagttgaaatgcactttcaagtaactcagatgtcactgaca
aacttactatagactggacatatgcccctccagcagcagccacacagtatcaatatttcattatcagcttctccagttacccaaccacagcaggcacatttgggatcgggatttccgtgggttgaaatgtatacaaaagggatgcattataagtaagaacccctaccataaaggacaatggg
acattcagctgtgctgtgaagaatccccagatgtgcatcataatattccatgacagagctaaacagtcacagaaaggggtttggcaccatgttctctgtggccctcttccatcctgtcttctgtgcccctcagccgtggtgtgtgctctgtgctgtgagatggggaggagggctgctg
ggctgaagaagaggagcaggtctgctataaagtcatttattgaggttccgatgacactgatcaggagggaagaggcgtgtatggcgaggcttgtgtccgtgtgctgagtgctggattcagactatgaagagacatattga