



Reference: PTXBC113589

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TMEM180; C10orf77; bA18I14.8; transmembrane protein 180; major facilitator superfamily domain containing 13A		

For research use only.

GENE DESCRIPTION	transmembrane protein 180
GENE ID	79847
ORIGIN SPECIES	Human
ACCESSIONS	BC113589
VECTOR	pDONR223
NCBI SYMBOL	TMEM180
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]