

PRR10-proline rich 10 Gene

Reference: PTXBC113593

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRR10; NCRNA00288; KIRREL3 antisense RNA 3		
For research use only.			

GENE DESCRIPTION	proline rich 10
GENE ID	283165
ORIGIN SPECIES	Human
ACCESSIONS	BC113593
VECTOR	pDONR223
NCBI SYMBOL	PRR10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaagcgcggagagtccttagacctcggagtgccgtgcaagaggggagaggagagagacgcccccctgctggaagccttcgctccggggcgccgcggggggaaggggtctctggacagtcgggtggagcggcagccctacagagacagccgaaagccagcaa
cttgaaaaccgccgaattctgggtttctgcacatccaggtggtccaggtctctgtgcccaccgtgcaagccgggtggatgccagcaggtggcatcctctccagcgtttcttcgcccagcgtggagtcggaagccgaatccatcagtcagcccccctcccaagcctccggtgc
caagtgctgggtctgtgagccctggccccgccccacctccgccagcggggcgagccggagctgggtgacccagaccctcgccgagcggggagcctacagagagtcggccggccccaggtttctgctggtggtgccaagaagtgacaggaaggcaaggctccccgg
aaagcgtccaagccctgctctctcttccaagccttgctgctgttgcggggggtttctgagatctgggacgaagagagtgatgtgcaaggttctgggtcccgggtgcagagtgagggtacggctgtcgagtgctcgagcaggtggcgtttggcacactgcagacctcagc
tcaccgccatttctctga