

## ATXN7L3-ataxin 7-like 3 Gene

Reference: PTXBC113595

|                               |                                          |                     |                |
|-------------------------------|------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                      | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                     | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | ataxin-7-like protein 3; ataxin 7 like 3 |                     |                |
| <b>For research use only.</b> |                                          |                     |                |

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | ataxin 7-like 3                                                                                   |
| <b>GENE ID</b>          | 56970                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC113595                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | ATXN7L3                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgaaaatggaggaaatgtcttggcctggataacagcaaaatagaggccatcgctcaggagatatacgggacctgtcaggattctgttgggattctgttgaggtaacacgggctgtcaagtgtggctacttcttggacgacacggaccctgatagcatgaaggatttga  
gatcgtggaccagccggcttggacatcttggacaggtttcaaccagtggagagcaaggagtggtttgccccaatgagtcgagcattgccctcccgcttggctcccatctggagaagtcctgggaatgggtcggacagcagccgaatcgaacccgggattgccaat  
agcaacaataatgaatagctgagagtgaccaagaagataatgatgacatcaatgacaacgactggctctatggctcggagaagaagccaagaagagaaagtcagacaagaaccccaattccctcgaagatccaagtcattaaacacaaaaatggggaacttagcaattcgatcct  
tttaagtatacaattcaactgggaltcagctatgagaccctggggccggaggagcttcgcagcctgtaaccacgcaatgtgggggtgatttctgaacacaccaagaagatgtgcacaaggtccctgcctgcccacagcacacagatgagcagaggcgaaccgtacggatttattctc  
ggccctcggctgtcctccagaggtcagagctccctggataatgacagcttggacatgactgacagccaggccctgatcagccggctcagtgaggacggctcctctgacctcaccctctgattcagctcctcaagacgagtgaatcagggatgggtctaggtaccaacagctct  
gagtcacggaaaaccaagaaaaagaatcccatctgagcctggttagggactgcctccggcctaggttccaacaagaagaagccaagccaccgacaccccgacgccagcatctatgatgacatcaactga