

PRR10-proline rich 10 Gene

Reference: PTXBC113597

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRR10; NCRNA00288; KIRREL3 antisense RNA 3		
For research use only.			

GENE DESCRIPTION	proline rich 10
GENE ID	283165
ORIGIN SPECIES	Human
ACCESSIONS	BC113597
VECTOR	pDONR223
NCBI SYMBOL	PRR10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaagcgcggagagtccttagacctcgagtgctgcaagaggggagaggagagagacgccccctgctggaagccttcgctccggggcgccgcgcgggggaagggtctctggacagtcggtggagcggcagccctacagagacagccgaaagccagcaa
cttggaagaccgcccgaattctgggtttctgcacatccaggctggctccaggctctgtgcccaccgtgcaagccggggatccagcagggtggcatctctccagcgtttcttcgccaggcgtggagtcggaggccgaatccatcagtcgccagccccctcccaagcctccggtg
ccaagtctgggtctgtgagccctggccccccccccacctccgccagcggggcgagccggagctgggtgacccagaccctcgcgaggcggggagcctacagaggtatccggccggccccaggttctgtctggtggccaagtgacagggaggcgaagctccccg
gaaagcgtccaagccctgctctctgctctccaagcctgctgctgtgtgccccagggtttctgagatctgggacgaagagagtgatgtgcaaggttctgggtcccggtgcaggagtgagggtacggctgctgagtgctcggagcaggctggcgtttggcacactgcagacctcag
ctcaccgccactttctctga