

CXorf1-chromosome X open reading frame 1 Gene

Reference: PTXBC113600

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CXorf1; transmembrane protein 257		

For research use only.

GENE DESCRIPTION	chromosome X open reading frame 1
GENE ID	9142
ORIGIN SPECIES	Human
ACCESSIONS	BC113600
VECTOR	pDONR223
NCBI SYMBOL	CXorf1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtattccagactattttatttaaaatcatcatatatactttgaacctttattttctaatgcaatcataaacttttaagtttttaataactagcatccccactgactatattctgttttgctttatctgctcaagcactttcgaccatattttatttaggattttttttatcttccacagctggatattgc
tatttcactttttatttcacatgcagctttaagacgtatgagcatcagcagcaaaatggttcagcctacagaatgcagctctccaggccctgccgaagaacatattctgtatgtctggccctacaaatga