



Reference: PTXBC113607

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDH20; FIB2; PCDH-GAMMA-B4; protocadherin gamma-B4; cadherin 20; fibroblast cadherin FIB2; fibroblast cadherin-2; protocadherin gamma subfamily B, 4		

For research use only.

GENE DESCRIPTION	protocadherin gamma subfamily B, 4
GENE ID	8641
ORIGIN SPECIES	Human
ACCESSIONS	BC113607
VECTOR	pDONR223
NCBI SYMBOL	PCDHGB4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr