



Reference: PTXBC113615

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GLUR2; GPRC1B; MGLUR2; mGlu2; metabotropic glutamate receptor 2; glutamate receptor homolog; glutamate receptor, metabotropic 2; glutamate metabotropic receptor 2		

For research use only.

GENE DESCRIPTION	glutamate receptor, metabotropic 2
GENE ID	2912
ORIGIN SPECIES	Human
ACCESSIONS	BC113615
VECTOR	pDONR223
NCBI SYMBOL	GRM2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]