

GRIN2B-glutamate receptor, ionotropic, N-methyl D-aspartate 2B Gene

Reference: PTXBC113618

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EIEE27; GluN2B; MRD6; NMDAR2B; hNR3; glutamate receptor ionotropic, NMDA 2B; N-methyl D-aspartate receptor subtype 2B; N-methyl-D-aspartate receptor subunit 3; NR3; glutamate [NMDA] receptor subunit epsilon-2; glutamate receptor subunit epsilon-2; glutama		

For research use only.

GENE DESCRIPTION	glutamate receptor, ionotropic, N-methyl D-aspartate 2B
GENE ID	2904
ORIGIN SPECIES	Human
ACCESSIONS	BC113618
VECTOR	pDONR223
NCBI SYMBOL	GRIN2B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr