



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RP1-21O18.1-kazrin Gene

Reference: PTXBC113621

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAZ; kazrin; kazrin, periplakin interacting protein		
For research use only.			

GENE DESCRIPTION	kazrin
GENE ID	23254
ORIGIN SPECIES	Human
ACCESSIONS	BC113621
VECTOR	pDONR223
NCBI SYMBOL	RP1-21O18.1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgaattccctgatgaagcgcgagacctgcgagaagcactgaacgtgtccaaagaattgccaccaggctcagcatctctctgggagtcgctgtctgtaccaagtgaaactcagcaggaggagccctccaggagcgcggcgccctctcgcgagacgcagaacattgacccctggtgtg
accaacccagcgggtctcgaatgtggcttcgagatcagactgacgtgaaggagtagcgcagacaacctgcacaacagccgcctccattggtctgtctgtctgtcggagcccaattcaatgccgaggccatggccactgccttcgcacccagggaagcacatctctccgcgagacactg
gcagagagatgagcgcgcctctcccccacacacccacacccctccggaggctgagctgttttggaaacgcgcctccagcgcgcctccagctcagcgcgcgagcgaaggaacacagcagcgcctctcaagtacaaaggctgcggcctgcctccctggggaagataggaa
gggcttcacacccaagaatgccgattctcatgactgactgtctcttcaaacgaagattgcgcgagagatgcgaccccccagcagcgcctggaacatgctccgtctggaagcgtlaaacacacgcctgaggtctcagctctgataa

More info: Please contact us for any question: contact@proteogenix.fr