

TATDN3-TatD DNase domain containing 3 Gene

Reference: PTXBC113638

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TatD DNase domain containing 3		

For research use only.

GENE DESCRIPTION	TatD DNase domain containing 3
GENE ID	128387
ORIGIN SPECIES	Human
ACCESSIONS	BC113638
VECTOR	pDONR223
NCBI SYMBOL	TATDN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgagcggctggcgttaggttggtggactgtcactgccactctccgccccggactttgacgcgatttggatgatgtgtggagaaagccaagaaggccaatgttggcccttggcagttgccgaacattcaggagaattgaaaagattatgcaacttcagaaaggtataatggg
tttgcctgccatgcttgggtgttcatccagttcaaggacttccaccagaagaccaaagaagtgtcacactaaaggatttggatgtagctttgccattattgagaattataaggatcgggtgttggcaattggagaggttggactagatttctccccagatttgcgcactggtaacagaagg
aagagcaaaagacaagtcctaatcagacagatccagtttagccaaaagactaaatttgcctgtaaatgtgcactcagctctgtggaagacctaaccatcaacctttacaagagcaaggtgctgagaaggtactgctgcatgcatgtggtcggccatctgtagccatggaaggagtaagag
ctgggtacttctcctcaattcccccttctatcataagaagtgacagaagcagaaccttggaaacaattgaccttaacttctatagcttagaacaagattcactgcactaggaccagaaaaacagggtacggaatgagccctggaacatttctattcagcagaatattgcccaggtgaaagg
gatctcagtggaagaagtatatagaagtgacacacagaatgcattaaaactgttcttaagctccgacactgtctccagaaatag