

IFNA4-interferon, alpha 4 Gene

Reference: PTXBC113642

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IFN-alpha4a; INFA4; interferon alpha-4; IFN-alpha-4; interferon alpha-4B; interferon alpha-76; interferon alpha-M1; interferon alpha 4		

For research use only.

GENE DESCRIPTION	interferon, alpha 4
GENE ID	3441
ORIGIN SPECIES	Human
ACCESSIONS	BC113642
VECTOR	pDONR223
NCBI SYMBOL	IFNA4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctgtccttttactgatggccgtgctggtcctcagctacaaatccatctgttctctgggctgtgatctgcctcagaccacagcctgggtaataggaggccttgatactctggcacaatgggaagaatctctatttctcctgcctgaaggacagacatgattcggattcccg
aggaggagtttgatggccaccagttccagaaggctcaagccatctctgtcctccatgagatgatccagcagaccttcaatctctcagcacagagagctcatctgtctgtgggaacagagcctcctagaaaaatttccactgaactttaccagcaactgaatgacctggaagcatgtgtgata
caggaggttggggtggaagagactccccgatgaatgaggactccatctggtctgtgaggaaacttccaaagaatcactctttatctaacagagagaagaatacagcccttgtgcctgggaggtgtcagagcagaatcatgagatccctctctgtttcaacaactgcacaaaaagattaa
ggaggaaggattga