

COL24A1-collagen, type XXIV, alpha 1 Gene

Reference: PTXBC113654

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	collagen alpha-1(XXIV) chain; collagen, type XXIV, alpha 1; collagen type XXIV alpha 1 chain		

For research use only.

GENE DESCRIPTION	collagen, type XXIV, alpha 1
GENE ID	255631
ORIGIN SPECIES	Human
ACCESSIONS	BC113654
VECTOR	pDONR223
NCBI SYMBOL	COL24A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr