



Reference: PTXBC113657

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GBD4; STSL; ATP-binding cassette sub-family G member 8; ATP-binding cassette, sub-family G (WHITE), member 8; sterolin 2; ATP binding cassette subfamily G member 8		

For research use only.

GENE DESCRIPTION	ATP-binding cassette, sub-family G (WHITE), member 8
GENE ID	64241
ORIGIN SPECIES	Human
ACCESSIONS	BC113657
VECTOR	pDONR223
NCBI SYMBOL	ABCG8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]